



| | |
|-----------------|---|
| test parameters | column: Ascentis Express Peptide ES-C18, 10 cm x 3 mm I.D., 2.7 µm particles (53313-U) |
| | mobile phase: (A) 0.1% (v/v) formic acid, pH 4.0 (titrated with ammonium hydroxide); (B) 50:50, (0.175% formic acid in water, pH 4.0) : acetonitrile |
| | gradient: 20 to 50% B in 15 min |
| | flow rate: 0.6 mL/min |
| | pressure: 207 bar |
| | column temp.: 35 °C detector: UV, 215 nm |