


1. Nivalenol
2. Deoxynivalenol
3. 15-Acetyldeoxynivalenol
4. 3-Acetyldeoxynivalenol



test parameters	column: Ascentis Express Phenyl-Hexyl, 10 cm x 2.1 mm I.D., 2.7 µm particles (53336-U)
	mobile phase: (A) water, acetonitrile, 90:10; (B) acetonitrile
	gradient: 0% B held for 0.5 min; to 15% B in 1 min; held at 15% B for 6 min
	flow rate: 0.3 mL/min
	pressure: 192 bar (2780 psi)
	column temp.: 30 °C
	detector: UV, 218 nm