


1. Valacyclovir (50 µg/mL)
2. Imiquimod (50 µg/mL)

Gradient Program		
Time (min)	%A	%B
0	95	5
10	20	80



test parameters	column: Ascentis C8, 15 cm x 2.1 mm I.D., 3 µm particles (581402-U)
	mobile phase: (A) 10 mM ammonium phosphate (pH 2.5 with phosphoric acid); (B) acetonitrile
	flow rate: 0.2 mL/min
	column temp.: 35 °C
	detector: UV, 254 nm
	injection: 2 µL
	sample: as indicated in mobile phase