


ニコチン酸誘導体

Nicotinic acid and derivatives

F000707B



Column : Hydrosphere C18 (5 μ m, 12nm)
150 $\times$ 4.6 mm I.D.
Eluent : 50mM KH₂PO₄ - K₂HPO₄ (pH 5.9) / acetonitrile (90/10)
Flow rate : 1.0 mL/min
Temperature : 37 $^{\circ}$ C
Detection : UV at 260 nm, 0.32 AUFS
Injection : 10 μ L (0.1 mg/mL)